

# LABOSHAKER R100

## Seesaw motion의 rocking mixer

- | Seesaw motion으로 부드럽고 균일한 mixing 구현
- | 다이얼 방식의 편리한 속도 조절
- | Set, Run, Remaining Time의 확인이 용이한 디지털 패널
- | 컴팩트한 구조로 최대 4개의 microwell plate 장착 가능
- | 별매의 엠보싱 매트와 고무 밴드로 다양한 용기를 안정적으로 사용 가능

## | Application

- Gel staining / destaining
- Western blot 반응 및 washing
- ELISA 반응을 비롯한 모든 Ag-Ab 반응
- 모든 종류의 mixing 및 반응 촉진에 유용





## Technical Specifications

Speed range (RPM)	0~60
Tilt angle ( ° )	±11
Max. load (Kg)	2.0
Platform area (W x D, mm)	270 x 176
Mains voltage (V)	220
Mains frequency (Hz)	50/60
Connected load (VA)	30
Dimensions (W x D x H, mm)	300 x 250 x 154
Run time	99(h):59(m):59(s), continuous
Weight (Kg)	5.8

## Ordering Information

Cat. No.	Product	Specification
GS-R100	LaboShaker R100 including Basic Plain Mat	220V, 50/60Hz
GSA-R1-M1	Basic Plain Rubber Mat of R100	270 x 176 x 3t
GSA-R1-M2	Non-slip Spiked Rubber Mat of R100	270 x 176 x 10
GSA-R1-B	Rubber Band	Φ120 x 2.62t Φ158 x 2.62t
GSA-R1-P2	Additional Platform for R100 (incl. Stacking support and Basic Plain Rubber Mat)	294.5 x 294.5
GSA-R1-P6	Stacking Support for R100 (1 each)	