

ARA



ARA mini 6

Features

- One spin action simply by closing the lid
- Safety stop by opening the lid
- Selectable operation modes between short spin and continuous 15 minutes
- Interchangeable rotors for 6 micro tubes and 2 PCR strips
- Electro sensitive non-contacting lid sensor

Technical Specifications

Max. RPM	6,500 RPM
Max. RCF	2,234 x g (1.5/2.0 mL Tube) / 2,182 x g (PCR Strip)
Capacity	6 x 1.5/2.0 mL, 2 x 8-Tube PCR Strip
ACC Time to Max Speed	3 sec
DEC Time from Max Speed	4 sec
Operation Mode	Selectable, Short Spin (10 sec) or Continuous (15 min)
Dimension (W x D x H)	150 x 150 x 90 mm
Weight	650 g
Power Supply	AC 110~220V, 50/60 Hz, 3A
Cat no.	AC-mini6