

ARA



Laboratory Shaker **D-ST**

Features

- 30 x 30 cm platform with an "anti-slip" robust shaker
- Max. load 3 kg
- Adjustable up to 140 rpm with tilting up to 15°
- Supplied with extra set of silicon pillar
- Optional non-slip spiked mats for 50 mL and 15 mL tubes
- Optional rubber bands for fixing containers

Technical Specifications

Speed Range	0 ~ 60 rpm
Shaking Angle	± 11°
Max. Load	2.0 kg
Platform Area (W X D)	270 x 176 mm
Mains Voltage	220 V
Mains Frequency	50/60 Hz
Power Consumption	30 VA
Dimensions (W X D X H)	300 x 250 x 154 mm
Weight	Timed < 100 hrs or continuous
Power Supply (V/Hz)	5.8 kg
Cat no.	D-ST