

# ARA



## Vortex Mixer **D-VM1**

### Features

- Available manual working mode (Touch on Mode)
- Adjustable speed upto 4,500 rpm
- Various kinds of accessories
- High quality and safety design
- CE certificate

### Technical Specifications

|                        |                   |
|------------------------|-------------------|
| Max. Speed             | 4,500 rpm         |
| Operation Mode         | Touch On          |
| Head                   | Mixing Cup Head   |
| Dimensions (W x D x H) | 100 x 100 x 67 mm |
| Net. Weight            | 550 g             |
| Cat no.                | D-VM1             |

## Vortex Mixer **VM2**

### Features

- Available manual working mode (Touch on Mode)
- Available automatic working mode (Constant Mode)
- Adjustable speed upto 3,000 rpm
- Various kinds of accessories
- High quality and safety design
- CE certificate

### Technical Specifications

|                        |                                     |
|------------------------|-------------------------------------|
| Max. Speed             | 3,000 rpm                           |
| Operation Mode         | Touch on or Continuous              |
| Dimensions (W x D x H) | 130 x 155 x 155 mm                  |
| Net. Weight            | 3 kg                                |
| Power                  | Universal(90 to 230V, 50/60Hz), 15W |
| Cat no.                | D-VM2                               |