


D-RT6 Roller Mixer



Features

- Constant speed control proportional to the medium or capacity
- Resistant to corrosion in organic solvents by black SUS 304 body
- Stable weight and durability for size

Technical Specifications

| | |
|--|--------------------------------|
| Number of Tube Rollers | 6 rollers |
| Maximum Capacity | 10 x 50 mL Tubes |
| Speed | Variable, 3 to 80 rpm |
| Speed Accuracy | ±5 |
| Speed Steps | 1 rpm |
| Rolling Angle & Tilting Angle | Rolling : ±3°, Tilting : ±5° |
| Time Control | Timed < 100 hour or Continuous |
| Dimension (W x D x H) | 300 x 326.5 x 135 mm |
| Weight | 7.5kg |
| Power Supply | AC 100~240V, 50/60Hz |
| Power Requirement | 9W |