


D-R8 Roller Mixer



Features

- Constant speed control proportional to the medium or capacity
- Resistant to corrosion in organic solvents by black SUS 304 body
- Stable weight and durability for size

Technical Specifications

Number of Tube Rollers	8 rollers
Maximum Capacity	14 x 50 mL Tubes
Speed	Variable, 3 to 80 rpm
Speed Accuracy	±5
Speed Steps	1 rpm
Rolling Angle & Tilting Angle	Rolling : ±3°
Time Control	Timed < 100 hour or Continuous
Dimension (W x D x H)	334 x 420 x 103.5 mm
Weight	6.1kg
Power Supply	AC 100~240V, 50/60Hz
Power Requirement	9W