


D-SR Laboratory Shaker



SR1



SR2

Features

- Reciprocal type
- Brushless DC motor with low noise (30 dB)
- Strong beltless drum without belt
- Stable working on 10 rpm and 999 hr integration timer
- Trustworthy safety device
- Easy to use various size

Technical Specifications

	D-SR1	D-SR2	D-SR3
Max. Speed		250 rpm (10 rpm/min)	
Max. Capacity	100 mL x 16, 250 mL x 9, 500 mL x 9, 1,000 mL x 4, 2,000 mL x 2	100 mL x 36, 250 mL x 25, 500 mL x 16, 1,000 mL x 9, 2,000 mL x 4	100 mL x 60, 250 mL x 40 500 mL x 24, 1,000 mL x 15, 2,000 mL x 8
Display		Digital LED Display	
Integral Timer		999.9 hr	
Drive and Motor		Beltless Direct Drive, Brushless DC Motor	
Plate Size	350 x 320 mm	450 x 450 mm	760 x 470 mm
Dimension (W x D x H)	360 x 360 x 110 mm	470 x 490 x 115 mm	745 x 490 x 115 mm
Weight	20 kg	28 kg	40 kg
Power Supply		AC 100~240V, 50/60Hz	
Power Requirement		30VA	