



D-ST Laboratory Shaker



1-layer



2-layer

Features

- 27 x 18 cm platform with an “anti-slip” robust shaker
 - Max. load 2 kg
 - Adjustable up to 60 rpm
 - Time control of <100 hr or continuous
 - Digital display of SET, RUN and remaining TIME
 - Optional non-slip spiked mats for 50 mL and 15 mL tubes
 - Optional rubber bands for fixing containers.
- The digital control with the display of SET, RUN and remaining TIME enables exact time control.

Technical Specifications

Speed Range	0 ~ 60 rpm
Shaking Angle	± 11°
Max. Load	2.0 kg
Platform Area (W X D)	270 x 176 mm

Power Supply	AC 100~240V, 50/60Hz
Power Requirement	30AV
Dimensions (W X D X H)	300 x 250 x 154 mm
Weight	5.8 kg