



# D-SD Laboratory Shaker



1-layer



2-layer

## Features

- 30 x 30 cm platform with an "anti-slip" robust shaker
- Max. load 3 kg
- Adjustable up to 140 rpm with tilting up to 15°
- Supplied with extra set of silicon pillar
- Optional non-slip spiked mats for 50 mL and 15 mL tubes
- Optional rubber bands for fixing containers

## Technical Specifications

<b>Speed Range</b>	0 ~ 140 rpm
<b>Shaking Angle</b>	0 ~ 15°
<b>Max. Load</b>	3.0 kg
<b>Platform Area (W X D)</b>	295 x 295
<b>Time Control</b>	99 (h) : 59 (m), continuous

<b>Power Supply</b>	AC 100~240V, 50/60Hz
<b>Power Requirement</b>	30 VA
<b>Dimensions (W X D X H)</b>	360 x 360 x 257
<b>Weight</b>	8.2 kg